
1
10
36b
36a
35b
35a
54
50
51
55
52
11
34b
34a
33b
33a
32b
32a
31b
31a
36
35
34
33
32
31
12
100
120
602
46
118
112
5

FIG. 1


100
102a
102b
102c
102d
140
111c
132
111b
180
111e
190
111f
178
111d
120
108
104
168
105
118
106
111a
110
112

FIG. 2A



FIG.2B


Bubble Trap
Cardioplegia Heat Exchanger
Bubble Sensor
Sequestration Reservoir
Venous Entry Module
PRE-BYPASS FILTER
Transfer Bag
Venous Reservoir
Arterial filter
Bubble Sensor
Oxygenator
Heat Exchanger
P1
P2
P3
P4
P5
P
120
301
106
122
110
112
114
116
118
119a
1196
122
126
128
129
132
134
136
138
140
142
144
146
148
150
152
153
156
158
160
162
164
166
168
170
172
174
175
176
178
180
186
188
190
192
192a
192b
302
303
304
305
306
307
308
309a
309b
310
311
320
322
330
332
334
336
338
340
342
344
346
348
350
352
354
356
358
360
362
364
366
368
370
372
374
376
378
380
382
384
386
388
390
392
394
396
398
400
401
402
403
404
405
406
407
408
410
412
414
416
418
420
422
424
426
428
430
432
434
436
438
440
442
444
446
448
450
452
454
456
458
460
462
464
466
468
470
472
474
476
478
480
482
484
486
488
490
492
494
496
498
500
502
504
506
508
510
512
514
516
518
520
522
524
526
528
530
532
534
536
538
540
542
544
546
548
550
552
554
556
558
560
562
564
566
568
570
572
574
576
578
580
582
584
586
588
590
592
594
596
598
600
602
604
606
608
610
612
614
616
618
620
622
624
626
628
630
632
634
636
638
640
642
644
646
648
650
652
654
656
658
660
662
664
666
668
670
672
674
676
678
680
682
684
686
688
690
692
694
696
698
700
702
704
706
708
710
712
714
716
718
720
722
724
726
728
730
732
734
736
738
740
742
744
746
748
750
752
754
756
758
760
762
764
766
768
770
772
774
776
778
780
782
784
786
788
790
792
794
796
798
800
802
804
806
808
810
812
814
816
818
820
822
824
826
828
830
832
834
836
838
840
842
844
846
848
850
852
854
856
858
860
862
864
866
868
870
872
874
876
878
880
882
884
886
888
890
892
894
896
898
900
902
904
906
908
910
912
914
916
918
920
922
924
926
928
930
932
934
936
938
940
942
944
946
948
950
952
954
956
958
960
962
964
966
968
970
972
974
976
978
980
982
984
986
988
990
992
994
996
998
1000

FIG. 3A


162
98
160
136
99
133
138
140
154
HYDROPHOBIC MEMBRANE
164
142
152
Bubble Trap
150
36
158
96
BUBBLE SENSOR
153
T3
Cardioplegia Heat Exchanger
Crystalloid Filter
16
93
143
144
142
18
148
146
38
134a
132
134
HEMO-CONCEN TRATOR
35
39
42
P3
172
176
P9
86
130
184
180
34
Waste Bag
Transfer Bag
37
44
P4
186
188
134b
190
33
194
14
P1
40
P2
170
179
46
VENOUS ENTRY MODULE
104
108
95
166
32
182
178
PRE-BYPASS FILTER
168
119
P7
175
93
VENOUS RESERVOIR
106
94
110
92
88
126
125
118
116a
Bubble Sensor
Arterial filter
Bubble Sensor
Oxygenator
P8
31
84
122
1166
91
90
116
114
112

FIG. 3B

TO9850 E62E9660



FIG. 4


108
530
544
536
538
540
542
546
534
532

FIG. 5A


554
530
536
108
544
546
556
534
532

FIG. 5B


541
543
539
548
550
83
552

FIG. 5D


85
560
548a
548b
548
558
550
81
552
552a
83
552b

FIG. 5C



FIG. 5E



FIG. 5F



FIG. 6A



FIG. 6B



FIG. 6C



FIG. 6D



FIG. 6E



FIG. 7A



FIG. 7B


500
105
506
505a
510
507
507a
525
508a
507b
512
504b
504a
502
502a
502b
112
501
513
508b
511
516

FIG. 8A



FIG. 8B



FIG. 8C

D 595 896 S 1 09/29/2009

FIG. 8D



FIG. 8E


522a
510
521a
508a
509a
520a
512
514
520b
521b
522b
517
516
511
520c
513
508b
522c
521c



FIG. 9A



FIG. 9B



FIG. 9C



FIG. 9D



FIG. 9E



FIG. 9F



FIG. 10A



FIG. 10B



FIG. 10C



FIG. 10D


728
727
Electronic Vacuum Regulator
721
726
Hydrophobic Filter
725
P
721
724
723
722
Venous Reservoir
106
111a
110
31
112

FIG. 11

D0959676 - 09-26-96



FIG. 12


120
150
176
188
174
140
132
301
180
190
311
303
129
310
119b
128
119a
178

FIG. 13


148
302
344
120
140
150
176
188
174
132
301
180
190
311
178
303
129
310
119b
128
119a

FIG. 14



FIG. 15



FIG. 16



FIG. 17



FIG. 18


B
149a
149b
404a
18a
89a
42a
40a
A
44a
A
403a
402a
405a
401a
14a
119a
119b
406a
129
B

FIG. 19A

Footage By Design



FIG. 19B



FIG. 20A



FIG. 20B



FIG. 21


149a
148
149b
404
307
140
89
18
42
176
40
142
301
790
132
280
282
188
164
185
180
44
190
272
176
402
182
403
130
405
401
309b
14
178
119a
128
119b
406
129

FIG. 22



Fig 23 A



Fig. 23B



FIG 24A



FIG 24B

TD9260 "E6/E9660



FIG. 25


340
358
Vaporizer (Optional)
350
348c
354c
344c
5%CO2 95%O2 Flow Meter
Pressure Regulator
Filter
5%CO2 95%O2 Source
352c
346c
342c
360
362
Filter
356
Total Flow Flow Meter
348b
344b
02 Flow Meter
Pressure Regulator
Filter
O2 Source
352b
354b
346b
342b
344a
342b
112
366
Liquid Leak Detector
Oxygenator
Pressure Regulator
Filter
Air Source
Scavange Line
364
352a
348a
346a
342a
Dryer Cross Flow
Sampling Pump
Sampling Pump
368a
368b
Flow Restrictor
Dryer
367
10
O2/CO2 Monitor
370
Optional Gas Concentration Monitoring System

FIG. 26


200
VLC Full Open
VLC 100 %
VLC Full Close
Venous
0
SAT XX%
HCT XX%
Temp. 0.0°C
Arterial
Line Pressure 0 mmHg
Flow 0.00 lpm
Temp. 0.0°C
Cardioplegia
Line Pressure 0 mmHg (Antegrade)
Coronary Sinus 0 mmHg
Flow 0.000 lpm
Bolus 0 ml
Total 0 ml
Temp. 0.0°C
Ischemic Time 00:00
Blender
Flow 3.2 lpm
FiO2 50%
FiCO2 0%
Other
Patient Arterial 110 mmHg
Patient Temp. 32.0°C
FeCO2 1%
8:47
220
242
User Set Up
Follow instructions and then press 'Load' to go to Load screen.
244
A-V
246
CPG
248
Suction/Fluids
250
Gases
252
Waveforms
254
Settings
1. Insert oxygenator and venous reservoir in holders.
2. Snap in pre-bypass filter and venous entry module in holders and place venous line in clamp and close cover.
3. Insert cartridge and arterial filter in holders.
4. Connect lines to A) arterial filter, B) venous entry module, and C) venous reservoir(2).
5. Place line in bubble sensor and close cover.
6. Place arterial and cardioplegia table lines in clamps and close covers.
7. Install pump loops and close all lids. Hang table pack on console.
Unload
Load
240
243

FIG. 27



FIG. 28A



FIG. 28B



FIG. 28C



FIG. 28D



FIG. 28E



FIG. 28F


222
224
226
228
230
220
VLC Full Open
222a
VLC 0%
222b
VLC Full Close
222c
Venous
1307
SAT
85%
HCT
26%
Temp.
21.3°C
Arterial
Line Pressure
151 mmHg
Flow
3.81 lpm
Temp.
21.0°C
Cardioplegia
Line Pressure
69 mmHg (Antegrade)
Coronary Sinus
0 mmHg
Flow
0.276 lpm
Bolus
360 ml
Total
360 ml
Temp.
20.9°C
Ischemic Time
00:00
Blender
Flow
3.2 lpm
FiO2
50%
FiCO2
0%
Other
Patient Arterial
110 mmHg
Patient Temp.
32.0°C
FeCO2
1%
8:55

FIG. 2B


242
User Set Up
A-V
CPG
Suction/Fluids
Gases
Waveforms
Settings
244
246
248
250
252
254
Follow instructions and then press 'Load' to go to Load screen.
245
1.
Insert oxygenator and venous reservoir in holders.
2.
Snap in pre-bypass filter and venous entry module in holders and place venous line in clamp and close cover.
3.
Insert cartridge and arterial filter in holders.
4.
A
B
C
Connect lines to A) arterial filter, B) venous entry module, and C) venous reservoir (2).
5.
Place line in bubble sensor and close cover.
6.
Place arterial and cardioplegia table lines in clamps and close covers.
7.
Install pump loops and close all lids. Hang table pack on console.
Unload
Load
256a
243

FIG. 30A



FIG. 30B

243



FIG 30C



FIG. 30D



FIG 30E



FIG. 30F



FIG. 30G



FIG. 30H



FIG. 30I



FIG. 30J



FIG. 30K



FIG. 30L



FIG 31A

FIG. 31B



FIG. 31C



FIG. 31D



FIG. 31E



FIG. 31F



246

246

FIG. 32A

FIG. 32B



FIG. 32C

FIG. 32D


267
268
Main
A-V
CPG
Suction/Fluids
Gases
Waveforms
Settings
Sensors
Zoom In
Zoom Out
Backward
Forward
252
243

FIG. 32E


Main
A-V
CPG
Suction/Fluids
Gases
Waveforms
Settings
254
Protocol
Sensors
CPG
Timers
Pulse
Other
284
284a
284b
284c
284d
284e
284f
286a
Current Protocol
Factory Default
286b
286b'
Wake-Up Protocol
Factory Default
286c
286c'
Load Protocol
Mike's Protocol
286
286d
286d'
Display Settings
Factory Default
286e
286e'
Save Protocol
Bill's Protocol
243

FIG. 33A



FIG. 33B



FIG. 33C



FIG. 33D



FIG. 33E


Main
A-V
CPG
Suction/Fluids
Gases
Waveforms
Settings
Protocol
Sensors
CPG
Timers
Pulse
Other
Bag Low Warning Alarm
ON
OFF
Bag Low Audio Alarm
ON
OFF
Bag Empty Pump Stop
ON
OFF
K+ High Warning Alarm
ON
OFF
K+ High Audio Alarm
ON
OFF
Select & Modify Bag Presets
Bag 1
Bag 2
Preset 1
500 ml 100 mEq/L
Modify
Preset 2
500 ml 100 mEq/L
Modify
Preset 3
500 ml 100 mEq/L
Modify
Preset 4
500 ml 100 mEq/L
Modify
Configure Bolus
Mode
Volume
500
Bolus Count
Up
Down
Select Delivery
Select Bag From Which To Deliver
1
2
Antegrade
Retrograde
Crystalloid only
Blood/Cryst. ratioed

FIG. 33F